


Vascular Intervention // Peripheral
Stainless Steel Guide Wire / 0.018"

Cruiser-18



Single-piece stainless steel core wire



Re-shapeable tip design with highly visible coil material



Multiple tip configurations offer options in complex lesions



Cruiser-18

Indicated for use in the peripheral vasculature.*

Technical Data		Guide wire		
		Diameter	0.018"	
		Length	195 cm and 300 cm	
		Material	Core wire Coil outer Coil inner	
			Stainless steel 6 cm stainless steel 2 cm platinum/tungsten (radiopaque)	
		Radiopaque tip length	2 cm	
		Coating	Proximal Distal	
			PTFE 13.5 cm hydrophilic	
		Range of tip flexibilities	Medium and stiff	
Ordering Information		Catalogue number	Support Level	Length
		357259	Medium	195 cm
		357260	Stiff	195 cm
		357261	Medium	300 cm
		357262	Stiff	300 cm

*Indication as per IFU

Manufacturer:
BRIVANT LTD.
Parkmore West
Business Park,
Galway, Ireland
Tel. +353 91 385 037
Fax. +353 9176 6598

Distributor:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com