

CRM

Cardiac Rhythm Management

Cardiac Resynchronization Therapy

Data Sheet

Selectra

Telescopic coronary sinus lead delivery system

- Integrated slittable serial valve
- Ergonomic handle design
- Optimal torque transmission
- Hydrophilic inner coating
- Reliable shape retention



Selectra

Technical data

Outer catheter	
Working length	45 cm; 55 cm
Inner diameter	2.44 mm (7.3 F)
Outer diameter	2.91 mm (8.7 F)
Material	PEBAX
Inner coating	Hydrophilic coating


Dilator	
Working length	53.5 cm; 63 cm

Inner catheter	
Working length	65 cm; 75 cm
Inner diameter	1.83 mm (5.4 F)
Outer diameter	2.31 mm (6.9 F)
Material	PEBAX
Inner coating	Hydrophilic coating
Tip angle	50°; 90°
Guide wire compatibility	0.014"-0.035"
Lead compatibility	< 5 F
Proximal connection	Female luer-lock

Ordering information accessory kit	
Selectra "Accessory Kit"	375518

Package contents:

- 12 cm³ syringe with luer-lock
- 4 transvalvular insertion tools
- 1 Seldinger guide wire 150 cm
- 1 torquer
- 1 Selectra splitter tool
- 2 bi-directional check-valves
- 2 sealing caps
- 2 1-way stop-cocks



Ordering information outer catheters	
Selectra "BIO2" 45 cm	375523
Selectra "BIO2" 55 cm	375524

Selectra "Amplatz 6.0" 45 cm	375519
Selectra "Amplatz 6.0" 55 cm	375520

Selectra "Extended Hook" 45 cm	375527
Selectra "Extended Hook" 55 cm	375528

Selectra "Hook" 45 cm	375529
Selectra "Hook" 55 cm	375530

Selectra "Multipurpose EP (MPEP)" 45 cm	375531
Selectra "Multipurpose EP (MPEP)" 55 cm	375532

Selectra "Multipurpose Hook (MPH)" 45 cm	375533
Selectra "Multipurpose Hook (MPH)" 55 cm	375534

Selectra "Right" 45 cm	375535
Selectra "Right" 55 cm	375536

Selectra "Straight" 45 cm	375537
Selectra "Straight" 55 cm	375521

Package contents:

- 1 outer guiding catheter
- 1 dilator for guiding catheter

Ordering information inner catheters	
Selectra "IC-50" 65 cm	375545
Selectra "IC-50" 75 cm	375546

Selectra "IC-90" 65 cm	375547
Selectra "IC-90" 75 cm	375548

Package contents:

- 1 inner catheter

