

Vascular Intervention // Coronary

# Galeo Pro

Coronary Guide Wire

Indicated for PTCA applications



- High tensile stainless steel core for excellent shape retention
- Enhanced torque response for unique steerability
- Hydrophilic coating to minimize friction
- Broad selection of tip stiffness and support levels

# Galeo Pro

## Coronary Guide Wire

Technical Data	Stent	
Diameter	0.014"	
Length	190 cm and 300 cm	
Core wire material	High tensile stainless steel	
Proximal coil	21.5 cm, stainless steel	
Distal coil	26 mm, palladium, radiopaque	
Proximal (shaft) coating	PTFE	
Distal coating	Hydrophilic (30 cm)	
Range of tip flexibilities	High Flexible (HF), Flexible (F), Medium (M)	
Range of support levels	Standard and ES (Extra Support)	
Shaft markers	92 cm and 102 cm from distal end	

Order Information	Catalogue number	Tip shape	Tip flexibility	Label color code	Support level	Length
	389781	Straight	High Flexible	■ Red	Standard	190 cm
	389782	J				
	389783	Straight	Flexible	■ Orange	Standard	190 cm
	389784	J				
	389785	Straight	Medium	■ Blue	Standard	190 cm
	389786	J				
	389787	Straight	Flexible	■ Yellow	Extra support	190 cm
	389788	J				
	389791	Straight	Flexible	■ Orange	Standard	300 cm

Galeo Pro is part of the BIOTRONIK coronary solutions portfolio, including:

- Stents: **Orsiro, PRO-Kinetic Energy, PK Papyrus** ■ Balloons: **Pantera Lux, Pantera LEO, Pantera, AngioSculpt**
- Guide Wires: **Cruiser, Magnum**

**For ordering please contact your  
local sales representative**

Manufacturer:  
Concert Medical, LLC  
77 Accord Park Dr.  
Norwell, MA 02061 -1605  
USA  
Tel +1 (781)871 7882  
info@concertmedical.com  
www.concertmedical.com

Distributor:  
BIOTRONIK AG  
Ackerstrasse 6  
8180 Bülach · Switzerland  
Tel +41 (0)44 8645111  
Fax +41 (0)44 8645005  
info.vi@biotronik.com  
www.biotronik.com