


Guide Wires

Angiography / Diagnostics

All Guide Wires are pre-coat PTFE and can be used for a variety of angiography procedures.




Code	Sizes	Tip Shape	Flexibility Area	Features
HR-3831	0.035"x150cm	J 3mm	Double Distal	Fixed Core, Finger Straightenable
HR-5521	0.035"x150cm	J 3mm	Single Distal	Movable Core
HR-3822	0.035"x150cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-3823	0.035"x150cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable, Firm Body
HR-3812	0.035"x150cm	Straight	Single Distal	Fixed Core
HR-3031	0.035"x180cm	J 3mm	Double Distal	Fixed Core, Finger Straightenable
HR-3922	0.035"x260cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-3912	0.035"x260cm	Straight	Single Distal	Fixed Core, Finger Straightenable
HR-4822	0.038"x150cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-4862	0.038"x150cm	Straight	Single Distal	Fixed Core
HR-4922	0.038"x260cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-4912	0.038"x260cm	Straight	Single Distal	Fixed Core

Qty/Box: 50

Hydrophilic

Hydrophilic Guide wires are designed to provide superior lubricity for enhanced control and steerability. Suited for a wide range of anatomical and clinical angiography procedures.



Code	Sizes	Guide Wire Stiffness	Tip Shape
HR-2851	0.018"x150cm	Standard	Angled
HR-2051	0.018"x180cm	Standard	Angled
HR-3811	0.035"x150cm	Standard	Straight
HR-3818	0.035"x150cm	Stiff	Straight
HR-5654	0.035"x150cm	Standard	Angled
HR-5655	0.035"x150cm	Stiff	Angled
HR-6511	0.035"x180cm	Standard	Angled
HR-6512	0.035"x180cm	Standard	Straight
HR-3631	0.035"x260cm	Stiff	Angled
HR-5641	0.035"x260cm	Standard	Angled
HR-5642	0.035"x260cm	Stiff	Straight
HR-5652	0.035"x260cm	Standard	Straight
HR-4851	0.038"x150cm	Standard	Angled

Flexible Tip Length: 3 cm

Qty/Box: 5