

PERIPHERAL INTERVENTION



# Hydrophilic Guidewire

## Micro Guide

# HYDROPHILIC GUIDEWIRE

## » Excellent Control

One-piece, super elastic nitinol core with distal flattened design to ensure the flexibility and controllability.

## » Durable and Lubricity

Covered by PU jacket and durable hydrophilic coating, to guarantee optimal pushability and trackability.

## » Safe and Clear

Coil ensures soft and atraumatic, distal 3cm radiopaque makes the wire highly visible under x-ray.

## SPECIFICATIONS

Guidewire O.D.	Tip O.D.	Length	Radiopaque Length	Core Material	Tip Shape
0.018"	0.018"	150cm	3cm, Pt. Coil	Nitinol	Straight Angled
0.021"		180cm			
0.025"					

## ORDERING INFORMATION

REF Code	Outer Diameter	Length (cm)	Curve	Tip Profile	Stiffness	Accessories
103A21815BY	0.018"	150	Angled	0.018"	Stiff	Shaping Needle Torquer
103A21818BY	0.018"	180				
103A22115BY	0.021"	150				
103A22118BY	0.021"	180				
103A22515QY	0.025"	150				
103A22518QY	0.025"	180				
103Z21815BY	0.018"	150	Straight			
103Z21818BY	0.018"	180				
103Z22115BY	0.021"	150				
103Z22118BY	0.021"	180				
103Z22515QY	0.025"	150				
103Z22518QY	0.025"	180				
103S21815BY	0.018"	150	45°			
103S21818BY	0.018"	180	70°			
103Q21815BY	0.018"	150				
103Q21818BY	0.018"	180				



@APT Medical

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. China

Email: international@aptmed.com Tel: +86 755-83591013

PN: ENG-Hydrophilic Guidewire-210X285X2P-20211015

