




Braidin™ MG

Steerable Guiding Sheath

ONE SHEATH, REACH ALL

Braidin™ MG steerable guiding sheath with steerable sheath tip offers precise maneuverability to help reach hard-to-access areas of the heart, which gives you a superior experience in catheter access and support.

- Bi-directional steering up to 180° makes it easier for access to complex anatomy
- Precise placement for different anatomies
- Braided structure provides superior curve retention even in a long time procedure
- Enhanced radiopacity for easier and more accurate location
- Ergonomic handle for easier operation



ORDERING INFORMATION

French Size	REF			Total / Working Length(cm)	Side Hole	Guidewire Information
	A Curve	B Curve	D Curve			
4F	990408201	990408202	990408203	96/82	--	0.035" × 180cm × J Tip
5F	990508201	990508202	990508203	96/82	--	0.038" × 180cm × J Tip
6F	990608201	990608202	990608203	96/82	--	0.038" × 180cm × J Tip
7F	990708201	990708202	990708203	96/82	--	0.038" × 180cm × J Tip
8F	990808201	990808202	990808203	96/82	--	0.038" × 180cm × J Tip
8F	990808221	990808222	990808223	96/82	2	0.032" × 180cm × J Tip
8F	990806821	990806822	990806823	82/68	2	0.032" × 180cm × J Tip
8.5F	990858201	990858202	990858203	96/82	--	0.038" × 180cm × J Tip
8.5F	990858221	990858222	990858223	96/82	2	0.032" × 180cm × J Tip
8.5F	990856821	990856822	990856823	82/68	2	0.032" × 180cm × J Tip
9F	990908201	990908202	990908203	96/82	--	0.038" × 180cm × J Tip
9F	990906821	990906822	990906823	82/68	2	0.032" × 180cm × J Tip
10F	991008201	991008202	991008203	97/82	--	0.038" × 180cm × J Tip
10F	991006721	991006722	991006723	82/67	2	0.032" × 180cm × J Tip
12F	--	991208202	991208203	97/82	--	0.038" × 180cm × J Tip
12F	991206721	991206722	991206723	82/67	2	0.032" × 180cm × J Tip

*Steerable Length: A Curve 26mm, B Curve 39mm, D Curve 52mm

Note: Only 68cm (working length) could be compatible with APT Medical 89cm transseptal needle, reference number: 94180890



@APT Medical

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. China

Email: international@aptmed.com Tel: +86 755-83591013

PN: ENG-Steerable Guiding Sheath-210X285X2P-20221214

