




Enlarger™


Inflation device

INFLATION DEVICE

- Ergonomic handle design for quick and easy inflation control
- Glow-in-the-dark pressure gauge for viewing in low-light environments
- 20ml syringe and 30atm max pressure fit all your dilatation procedures



Lock



Release

ORDERING INFORMATION

REF	Volume(ml)	Pressure(atm/psi)	Stopcock
191020S	10	20/294	•
191020	10	20/294	
192020S	20	20/294	•
191020	20	20/294	
192030S	20	30/440	•
192030	20	30/440	
193030S	30	30/440	•
193030	30	30/440	



@APT Medical

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. China

Email: international@aptmed.com Tel: +86 755-83591013

PN: ENG-Inflation Device-210X285X2P-20211015



www.apptmed.com