

PERIPHERAL INTERVENTION



**Snare Kit**

# ONE SNARE, FITS ALL

For retrieving venous filters and foreign objects in the vasculature

» **Excellent visibility and shape retention**

Loop: Gold plated tungsten coil over Nitinol wires.  
Delivery Catheter: Distal Pt. Marker for precise locating.

» **Minimize friction and maximize torqueability**

Pusher: PTFE over Nitinol shaft.



## ORDERING INFORMATION

REF.	Loop Diameter(mm)	Snare Length(cm)	Catheter O.D.(mm)	Catheter length(mm)
57021750*/57021751	2	175	1.06 (3F) /-	150 /-
57022000*/57022001	2	200	1.06 (3F) /-	150 /-
57041750*/57041751	4	175	1.06 (3F) /-	150 /-
57042000*/57042001	4	200	1.06 (3F) /-	150 /-
57051250*/57051251	5	125	1.40 (4F) /-	150 /-
57071750*/57071751	7	175	1.06 (3F) /-	150 /-
57072000*/57072001	7	200	1.06 (3F) /-	150 /-
57100650*/57100651	10	65	1.40 (4F) /-	150 /-
57101250*/57101251	10	125	1.40 (4F) /-	150 /-
57151250*/57151251	15	125	2.06 (6F) /-	150 /-
57201250*/57201251	20	125	2.06 (6F) /-	150 /-
57250650*/57250651	25	65	2.06 (6F) /-	150 /-
57251250*/57251251	25	125	2.06 (6F) /-	150 /-
57301250*/57301251	30	125	2.06 (6F) /-	150 /-
57351250*/57351251	35	125	2.06 (6F) /-	150 /-

\* Kit with delivery catheter



@APT Medical

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. China

Email: international@aptmed.com Tel: +86 755-83591013

PN: ENG-Snare Kit-210X285X2P-20211015

