




PERIPHERAL INTERVENTION



Susrail™
Guidewire

SUSRAIL™ GUIDEWIRE

» Optimal Support

Stainless steel core provides excellent support and torqueability for device delivery in interventional procedures, including endovascular aortic repair.

» Core-to-Tip structure

Shapeable tip with good maneuverability.

» Duo-coil Design


Atraumatic tip with superior tactile feedback.

» Precise Visibility


Platinum tungsten radiopaque allows clear visibility.

» Durable Lubricity

Distal hydrophilic coating and proximal PTFE coating for crossing lesions smoothly.



ORDERING INFORMATION

REF.	Diameter	Length (cm)	Tip Shape
42135261	0.035"	260	Straight(S)
42635261	0.035"	260	Double J 



@APT Medical 🔍

Address: Room 601, Building B, Tongfang Information Port, No.11 Langshan Road, Nanshan District, Shenzhen, P.R. China

Email: international@aptmed.com Tel: +86 755-83591013

PN: ENG-Guidewire-210X285X2P-20230113

 **APT Medical**
www.aptmed.com