



PRODUCT CATALOGUE

DEFIBRILLATION

+EURODEFIPADS®+



EURODEFIPADS®

MULTIFUNCTION PADS



Manufactured with high-quality materials, subjected to strict and continuous process control and thoroughly tested for compatibility and compliance with international standards, EURODEFIPADS® ensure both safety and excellence.

Easy to use, they deliver optimal clinical performance in both emergency situations and elective procedures.

EURODEFIPADS® is a complete line of adult and pediatric disposable multifunction electrodes for defibrillation, synchronized cardioversion, external cardiac stimulation and ECG monitoring with different connectors compatible with the main defibrillators present in the market.

Radiotransparent models (with radiolucent cable and pads) have been designed for use in EP labs, haemodynamics and during radiographic and angiographic procedures. Lead-out preconnected models (marked with the suffix W) have the cable and connector placed outside the pouch, so the pads can be connected to the defibrillator before opening the packaging. This design is especially useful in emergency situations.


EURODEFIPADS® are classified as Class IIb medical devices in accordance with the Medical Device Regulation (EU) 2017/745 (MDR) and are CE marked following assessment by a Notified Body, ensuring full compliance with the highest standards of safety, performance and quality.

ADULT

PAEDIATRIC

RADIOTRANSSPARENT

LEAD OUT



EURODEFIPADS®



Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

CONNECTORS

BRANDS

ADULT

UNIVERSAL

PAEDIATRIC

Radiopaque

Radiotransparent

Radiopaque

AED

Radiotransparent

		ADULT		UNIVERSAL	PAEDIATRIC		
		Radiopaque	Radiotransparent		Radiopaque	AED	Radiotransparent
F7950P F7950PW F7750P							
		Philips Medical®					
F7950PL F7950WPL others							
F7951 F7951W F7951/300							
		Zoll Medical Corp® M&B					
F7952PL F7952WPL others							
		Physiocontrol® Bexen Cardio Cardioline® Mindray Comen Jousing Acoresmed					
F7953 F7953W							
		Agilent Laerdal Medical® Hewlett-Packard					
F7954 F7954W							
		Schiller Esaote					
F7955 F7955W							
		Nihon Kohden Welch Allyn® Artema Innomed MRL Artema S&W Dräger					

EURODEFIPADS®

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

CONNECTORS

BRANDS

ADULT

UNIVERSAL

PAEDIATRIC

Radiopaque

Radiotransparent

Radiopaque

AED

Radiotransparent

 F7956 F7956W F7756	 F7956P F7956PW F7756P	Schiller	 F7956 F7956W	 F7756		 F7956P F7956PW		 F7756P
		Cardiac Science GE®	F7957W				F7957PL	
		CU Medical Systems, inc	F7958 F7958W			F7958P F7958PW	F7958PL	
		CU Medical Systems, inc Cmos Drake Amoul® Ambulanc (Shenzhen) TECH. CO., LTD	F7959 F7959W	F7759		F7959P	F7959PL F7959WPL	
		Metrax Primedic	F7960W	F7760W				
		Nihon Kohden		F7761	F7961U F7961UW F7961UWXL F7761U	F7961P F7961PW		
		Marquette GE® Esaote	F7962	F7762		F7962P		
		Nihon Kohden	F7963					
		Weinmann	F7964W			F7964PW		
		Corpuls	F7965			F7965P		
 F7965W/NG F7965/PLUSNG	 F7965PW/NG	Corpuls	F7965W/NG F7965/PLUSNG			F7965PW/NG		

EURODEFIPADS®

Please contact the International Sales Dept. in order to get the compatibility list and find out the right reference for each defibrillator brand and model.

Brands, models and commercial names reported in this document belong to the companies owning the corresponding registered brands, models and names. FIAB is not affiliated to any of the above mentioned companies.

CONNECTORS		BRANDS							
<div>F7966W</div> 	<div>F7966WPL</div> 	Defibtech	F7966W				F7966WPL		
<div>F7967 F7967W</div> 	<div>F7967P F7967PW</div> 	Schiller Philips	F7967 F7967W			F7967P F7967PW			
<div>F7969W</div> 	<div>F7969PW</div> 	Mediana	F7969W			F7969PW			
		GE®	F7981			F7981P			
		Mediana Tecnogaz	F7986 F7986W		F7986UW F7986UW/D700 F7786UW/D700	F7986P F7986PW			
<div>F7987W</div> 	<div>F7987PW</div> 	Zoll Medical Corp®	F7987W			F7987PW			
		Nihon Kohden Not for sale in France, Germany and UK			F7989UW				
		Radian			F7992U				
		Laerdal Medical®	F7970			F7970P			
		Physiocontrol®	F7971						
		Marquette GE® Esaote	F7972			F7972P			

Besides **EURODEFIPADS®** models, some additional models are intended for external cardiac stimulation only.



CONNECTORS	BRANDS		
	Zoll Medical Corp® (M series, pacer version only)	adult paediatric	F7951/STIM F7951P/STIM
	Physiocontrol®	adult	F7973
	Nihon Kohden Innomed	adult	F7955/STIM

TRAINER PADS

Defibrillation training pads (adult and paediatric) are equipped with snap connection and must be used with the appropriate connection cable. The pad's special shape makes the application extremely easy both on manikins and on other surfaces.

PADS

REF.	PATIENT
F7999	adult
F7999P	paediatric



CONNECTING CABLES

CONNECTORS	BRANDS	
	Philips Medical® Laerdal Medical® Agilent Philips Medical	F5350
	Physiocontrol® Osatu Bexen Cardio Cardioline® - Mindray	F5352
	Welch Allyn	F5355
	Schiller	F5356
	Weinmann Cardiaid	F5364
	Defibtech	F5366
	Mediana	F5369
	Mediana Tecnogaz	F5386
	Nihon Kohden	F5389



FIAB SpA
Via Costoli, 4 - I-50039 Vicchio (Florence) - Italy
Tel. +39 0558497999
export@fiab.it - www.fiab.it

