



3Flow

Aspiration Catheter

Indicated for central and peripheral circulatory system including saphenous vein grafts



- Extra large 60 ml locking syringe allows for longer aspiration time
- Large inner diameter enables extraction of larger thrombotic particles
- Hydrophilic distal coating for excellent trackability and crossability

3Flow

Aspiration Catheter

Technical Data	Aspiration catheter	3Flow 6F	3Flow 7F
	OD (dist./mid/prox.)	0.067"/0.067"/0.051"	0.078"/0.078"/0.063"
	Guide catheter min. ID	0.071" (1.80 mm)	0.081" (2.06 mm)
	Usable length	145 cm	145 cm
	Proximal shaft	PEEK	PEEK
	Distal shaft	SCP [Semi Crystalline Polymer]	SCP [Semi Crystalline Polymer]
	Coating (distal 25 cm)	Hydrophilic	Hydrophilic
	Distal marker band	Platinum-iridium, 3 mm from tip	Platinum-iridium, 3 mm from tip

Order Information	Catalogue number	Size
	387456	6F
	387457	7F

Manufacturer:
ARTHESYS
4, rue René Razel
91400 Saclay
France
www.arthesys.com

Distributor:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach · Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com

For ordering please contact your
local sales representative